



THE ASCLEPION EFFECT

MicroPulse Scanner



CO₂ fractional treatment -
as flexible as you need it

Fast, safe, re-producible, strong

SKIN REJUVENATION

REMOVAL OF PIGMENTS

FACIAL TIGHTENING

SMOOTH WRINKLES

REDUCTION OF PORE SIZES

IMPROVEMENT OF ACNE SCARS

EXCLUSIVE for MultiPulse

Technology made in Germany



Asclepion
Laser Technologies

MicroPulse Scanner

USPs/Benefits:

- Flexible cover rate
- Changeable spot pattern and sizes
- Clear protocols clinical proved
- FDA cleared
- No consumables



Before treatment



Before treatment



After treatment



After treatment

MultiPulse for a broad range of applications:

- Standard and special handpiece for dermatology, gynecology, ENT
- Scanners, Micromanipulators and Microscopes adapters
- Upgradeable

"As flexible as you need it" Adjustable parameters:

- Power (max. 30W)
- Fractional density (Pitch 200 - 2000 μ m)
- Dwell time (max. 2000 μ s)
- Ultrafast scan mode (normal, interlaced, autofill)
- Scan pattern (Shape, Ratio, Size)
- Adjustable repetition rate

Parameters specially adaptable to needs & wishes of patient



Display