

# SPECIFICATIONS\*

|                |                                             |
|----------------|---------------------------------------------|
| LASER          | MCL31 Dermablade (Er:YAG), class            |
| WAVELENGTH     | 2940 nm                                     |
| FLUENCE        | Max. 250J / cm <sup>2</sup>                 |
| PULSE DURATION | 100 – 1000 µs                               |
| FREQUENCY      | 1 – 20Hz                                    |
| MODES          | Ablation, Therm, MicroSpot                  |
| HANDPIECES     | 45° reflector, MicroSpot 9x9 mm             |
|                | 45° reflector, focus opt. 4.5 mm (optional) |

\*all specifications are subject to change without notice



**Asclepion**  
Laser Technologies



**MADE IN GERMANY**

Asclepion Laser Technologies GmbH • Brüsseler Str. 10 | 07747 Jena | Germany | [www.asclepion.com](http://www.asclepion.com)

MCL31 IV 001 0914 EN



# TURN BACK YOUR BIOLOGICAL CLOCK

WITH THE NEW  
**Juliet**  
*the feminine laser*



**INDICATIONS:**  
VAGINAL LAXITY  
STRESS INCONTINENCE

**Asclepion**  
Laser Technologies